



by David Glaser, Au.D., CCC-A, FAAA
Doctor of Audiology

GET WITH THE PROGRAM

As we go through our days at home, at work, and at play, we subject our ears to a number of different listening environments. We might have a quiet conversation over coffee in the morning, while work may be conducted in offices with multiple distractions and steady background noise. Because one program on a hearing instrument cannot be expected to perform equally well in different listening environments, many hearing instruments are outfitted with multiple programs. Each program may be accessed with the push of a button on the instrument or with the use of a remote control. More sophisticated hearing instruments are capable of analyzing the listening environment and adjusting automatically, without touching the hearing instrument or a remote control.

TIP OF THE WEEK

Some hearing instruments are so sophisticated that they can be adjusted directly through the cellular phone.

P.S. Some of the remote controls used to access programs on a hearing instrument can act as Bluetooth streaming devices

AA-48